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General directions


Infusion filter for Human Albumin 20% solution (HA); for use during infusion or to prepare extracorporeal HA treatment in liver disease in order to deplete conservatives (caprylate and/or acetyltryptophan).

CAVE! Use only with enclosed particle filter (RoweFil 48); suitable for use with infusion lines for infusion pumps or commercial infusion systems..

CAVE! Use at least 100 ml of isotonic sodium chloride per Hepalbin®-Adsorbent for pre rinsing!

Examples of application

Depending on the amount (quantity) of HA use of 6 Hepalbin®-Adsorbents in one application is possible. As in example B you can put 2 up to 4 Hepalbin®-Adsorbents in a row. If using 6 Hepalbin®-Adsorbents (example C) you need to add a Y-connector (two enclosed in a unit) in order to generate two columns. For use of 1 filter in 100 ml 20% albumin no infusion pump is required. If more than 1 filter is used in a row, an infusion pump must be used.

**Legend:**

- 1 Human Albumin 20% solution (HA) or sodium chloride isotonic 0.9% for pre-rinse and post-purge
- 2 infusion lines 2a without, 2b with infusion pump
- 3 Hepalbin®-Adsorbent
- 4 particle filter RoweFil 48
- 5 A, B venous access; C receptacle bag for extracorporeal treatment
- 6 Y-connector
- 7 3-way stopcock

Table: Recommendations with respect to the quantity of HA

HA (20%) volume [ml]	100	200	300	400	500	600	700	800
Number of adsorbents	1	2	3	3	4	4	6	6
Flow rate [ml/h]	150	300	450	450	600	600	900	900
Minimum of sodium chloride for post-purge [ml]	100	150	200	200	250	250	350	350

Preparation

- use a common metallic holder for infusion bottles for one or more Hepalbin®-Adsorbents
- assemble components as shown in figures (consider flow direction)
- avoiding bubbles, fill and rinse infusion system with sodium chloride (**CAVE!** Hepalbin®-Adsorbent and particle filter bottom up)
- during filling an rinsing, 100 ml sterile (infusion grade) 0.9% sodium chloride solution per adsorbent will enter the system out of which 30 ml to 50 ml will be discarded and 70 ml to 50 ml will remain in the system before connecting the albumin bottle
- discard excessive volume, attach HA, connect to patient

Performance

- **CAVE!** Close infusion line below particle filter in case of interruption such as changing bottles
- 1 drop per second (to be monitored in drip chamber) is the optimal flow rate without using infusion pump
- the minimum of a rinsing volume to wash out albumin remains out of the system is given in the table
- during filling and rinsing, 100 ml sterile (infusion grade) 0.9% sodium chloride solution per adsorbent will enter the system out of which 30 ml will be discarded if only 1 Hepalbin®-Adsorbent is used, but if more adsorbents are used, an additional volume of 50 ml per adsorbent shall be discarded; the rest of the rinsing solution remains in the device

Disposal

Discard Hepalbin®-Adsorbent complying legal requirements for hazardous waste.

Restrictions

Use only with Human Albumin 20% solution and infusion solutions which are approved as pharmaceuticals.

Precautions

- for single use only (the ability of binding toxins is limited and not reversable, reesterilisation leads to material damage!)
- pressure-resistant up to 1,2 bar
- sterile and free of pyrogens only if package is intact
- do not use after expiration date
- handle with hygienic operating principles

Performance characteristics

If applying as advised (table) there could be reached a ratio stabilizer/albumin well below 0.3 mol/mol.

Sterilization

The Hepalbin®-Adsorbent is sterilized by gamma radiation.

Storage

Storage of Hepalbin®-Adsorbent at room temperature, in dry conditions without direct sun exposure.