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This paper gives additional informations about practice and examples for application of the Hepalbin®-Adsorbent.

Table of contents

1 Purpose.....	1
1.1 Indication.....	1
1.2 Contraindication.....	2
2 Specification.....	2
3 Recommendations for various volumes of HA-solutions.....	2
3.1 Characteristics in extracorporeal therapy	2
4 Examples for applications.....	2
4.1 Example A, for the use of 1 Hepalbin®-Adsorbent for IV infusion	3
4.2 Example B, for the use of 2 to 4 Hepalbin®-Adsorbents in a row	3
4.3 Example C, for the use of 6 Hepalbin®-Adsorbents, e.g. for Single Pass Albumin Dialysis (SPAD).....	4
4.4 Example D, for the use of 3 Hepalbin®-Adsorbents in a row to purify 80g of Albumin for Extracorporeal Albumin Dialysis (ECAD).....	4
4.5 Technical Support Devices.....	5
5 Limitations.....	5
6 Warnings.....	5
7 Data.....	6
8 Troubleshooting.....	6

1 Purpose

Infusion filter for Human Albumin **20%** solution (**HA**); for use during infusion **or** to prepare extracorporeal HA treatment in liver disease in order to deplete stabilizers (caprylate [octanoate] and/or N-acetyltryptophanate).

1.1 Indication

HA is quantitatively the biggest part of the total protein in the plasma. When infused, the effect is to expand the circulating blood volume.

Its protective function is achieved by binding and releasing of endogenous and exogenous toxins, drugs, fatty acids, vitamins, hormones etc. Thus, HA takes control of free, hence active substances while simultaneously offering a reservoir.

In liver disease, HA is administered intravenously as a plasma expander and carrier for toxins or it is used in extracorporeal detoxification as a substitute in plasmapheresis or it is added to the dialysate solution in Albumin Dialysis in order to increase the removal of Albumin bound toxins.

However, as commercial Albumin undergoes various steps of virus inactivation and storage, industry is required to add stabilizer molecules to pharmaceutical preparations of Albumin. The purpose of said stabilizers is to bind as ligands at the benzodiazepine binding site of Albumin (also known as Sudlow II) and change the conformation to a more thermo-resistant protein molecule.

Unfortunately, this required process of stabilization for sterilization and storage has an undesired adverse effect on the theoretical benefit of providing additional binding sites to the patients, as the occupation of binding sites by said stabilizing ligands, which are caprylate (octanoate) and acetyltryptophanate, reduces the availability of binding sites.

In addition, research conducted after the process of Albumin stabilization had been established, revealed that tryptophan analogues and octanoate accumulate in liver disease and can precipitate vasodilation and hepatic coma in liver failure.

The purpose of the Hepalbin®-Adsorbent is to remove octanoate and N-acetyltryptophanate from Albumin when applied to patients with liver disease. This will increase the efficacy of Albumin therapy as the adverse effects of octanoate and tryptophan metabolites will not reduce the clinical benefit of Albumin.

The use of Albumin in the context of liver disease is established in the following indications:

- Paracentesis of Ascites
- Spontaneous Bacterial Peritonitis
- Hepatorenal Syndrom
- Albumin Dialysis or Apheresis (MARS, SPAD, etc.)

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1.2 Contraindication

No contraindications are known for using the Hepalbin®-Adsorbent. As its purpose is to be used in Albumin infusions, contraindications may arise from the contraindications of Albumin only.

2 Specification

Diameter:	108 mm
Height:	56 mm
Weight:	58 g
Material of shell	Polycarbonate
Fluid uptake:	50 ml
Durability:	59 months
Sterilisation:	Gamma radiation
Storage:	at room temperature, protect from light

3 Recommendations for various volumes of HA-solutions

One Hepalbin®-Adsorbent can significantly remove the stabilizers of 100 ml HA 20% with a flow rate of 150 ml/h. The efficacy of stabilizer depletion is reduced at higher flow rates but may still be clinically efficient when double flow rates are applied (**Fig. 5**).

The purification of more than 100 ml Albumin 20% requires the application of 2 to 6 Hepalbin®-Adsorbents and the use of an infusion pump.

Regarding the flow rate, one should abide by the rules according to Table 1. The flow rate of pre-rinsing can average at about 999 ml/h.

HA (20%) volume [ml]	100	200	300	400	500	600	700	800
Number of adsorbents	1	2	3	3	4	4	6	6
Flow rate [ml/h]	150	300	450	450	600	600	900	900
Minimum for pre-rinse [ml]	100	200	300	300	400	400	600	600
Minimum for post-rinse [ml]	100	150	200	200	250	250	350	350

Table 1 recommendation for depleting of stabilizers in different quantities

3.1 Characteristics in extracorporeal therapy

If the Hepalbin®-Adsorbent is used to prepare Albumin solutions for extracorporeal therapies in larger quantities (e.g. between 400 and 800 ml), it is frequently desired that the volume of the filtered product is not too much diluted or expanded by the pre- and post-rinsing procedures. One reason is that certain methods of Albumin Dialysis (MARS) require a concentration between 10 and 20% and another, that the volume of dialysate circuits may be limited (MARS Albumin circuit holds less than 800 ml).

In order to achieve this goal, specific standard operational procedures have been developed, which include a discarding of a defined amount of pre-rinsing solution while starting perfusion of Albumin and by recovering Albumin from the filter by defined post-rinsing procedures.

Using an infusion pump without a drop counting control mechanism may be preferred, as albumin tends to cause foam in the dripping chamber which may cause false alarms due to the missing signal induced by the drops.

The procedure will be described by giving the example D (4.4).

4 Examples for applications

Depending on the volume of HA 20% solution to be filtered (100-800 ml), up to 6 Hepalbin®-Adsorbents can be used in one application.

Using an infusion pump, 2 up to 4 Hepalbin®-Adsorbents can be applied in a row (see Fig. 1B or 1D). If 6 Hepalbin®-Adsorbents are required (for instance to purify 800 ml of albumin 20%, see table 1), two rows of 3 Hepalbin®-Adsorbents should be used in parallel, using a "Y-connector", in order to maintain the pressure within the acceptable range of a standard infusion pump.

For use of one adsorbent in 100 ml HA 20% no infusion pump is required. If using more than one adsorbent in a row, a infusion pump must be used.

Regarding the flow rate, one should abide by the rules according to Table 1. The flow rate of pre-rinsing can average at about 999 ml/h.


Figure 1: example of use A, B, C und D

Legend:

- 1: HA 20% solution or sodium chloride isotonic 0,9% for pre- and post-rinse
- 2: infusion lines 2a without, 2b with infusion pump
- 3: Hepalbin®-Adsorbent
- 4: particle filter RoweFil 48
- 5: A, B venous access; C, D collection bag for extracorporeal treatment
- 6: Y-connector
- 7: 3-way stopcock

- A: infusion of 100 ml HA (with or without infusion pump);
 B: infusion of 200 ml HA;
 C: preparation of 600 to 800 ml HA
 D: preparation of 300 to 400 ml HA

4.1 Example A, for the use of 1 Hepalbin®-Adsorbent for IV infusion

This example describes the use of 1 Hepalbin®-Adsorbent for the purification of 100 ml Albumin 20% solution and is pictured in **Fig. 1A**.

A Hepalbin®-Adsorbent (3), a particle filter (4) and an infusion line (2) are needed. The use of an infusion pump is recommended, but not required. As hospitals use different infusion systems, it is not part of the Hepalbin®-Adsorbent package.

1. Connect part 1 (appr. 100 ml NaCl), 2 (infusion line without or with an infusion pump), 3 (Hepalbin®-Adsorbent) and 4 (particle filter) as pictured in **Fig. 1A**.

2. During filling and pre-rinsing with 100 ml sodium chloride, the Hepalbin®-Adsorbent and the particle filter should be placed in a way that the fluid enters at the bottom and flows upwards. For achieving this, one could either use a plastic or metal clamp to hang the connected particle filter to the hook of the infusion stand which is usually situated at 1 meter height. This ensures a complete and almost bubble free filling of the Hepalbin®-Adsorbent and the particle filter. The fluid leaving the particle filter should be discarded until 100 ml are used. As the "dead volume" of one adsorbent and the infusion lines together account for approximately 70 ml, the volume of wasted sodium chloride during the rinsing should be around 30 ml if 100 ml sodium chloride are used for filling and rinsing.

3. Replace the sodium chloride with the 100 ml container of HA 20% solution. Connect the infusion line of the particle filter with the IV access of the patient 5 (e.g. a flexible needle or the connector for an IV catheter). Adjust the flow rate to approximately 150 ml/h (for use without infusion pump and infusion by gravity, the roller clamp of the infusion line should be adjusted that 1 drop per second is falling in the dripping chamber).

4. After the HA container is empty, replace it with 100 ml sodium chloride 0.9% and rinse the infusion system with up to 100 ml in order to recover the remaining albumin from the infusion system without changing the flow rate.

4.2 Example B, for the use of 2 to 4 Hepalbin®-Adsorbents in a row

This example describes the use of 2 Hepalbin®-Adsorbents for the purification of 200 ml Albumin 20% solution. The application is pictured in **Fig. 1B**.

Two Hepalbin®-Adsorbents (3) and one particle filter (4) are needed. Also, one infusion line and an infusion pump are required.

1. Connect part 1 (200 ml sodium chloride), 2 (infusion line with an infusion pump), 3 (2 Hepalbin®-Adsorbents) and 4 (particle filter) as pictured in **Fig. 1B**.

2. During filling and pre-rinsing with the sodium chloride amount specified in table 1, the Hepalbin®-Adsorbent and the particle filter should be placed in a way that the fluid enters at the bottom and flows upwards. For achieving this, one could either use a plastic or metal clamp to hang the connected particle filter to the hook of the infusion stand which is usually situated at 1 meter height. This ensures a complete and almost bubble free filling of the Hepalbin®-Adsorbent and the particle filter. The fluid leaving the particle filter should be discarded until 200 ml are used, according to table 1.

3. Replace the sodium chloride with 200 ml of HA 20% solution (If more than one container is needed, containers can be replaced consecutively.). Connect the infusion line of the particle filter with the IV access of the patient 5 (e.g. a flexible needle or the connector for an IV catheter). Adjust the flow rate to approximately 300 ml/h. If the remaining sodium chloride in the system should not be infused or used, up to 100 ml can be discarded before the

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infusion line should be connected to the patient. This will reduce the dilution without risking significant loss of Albumin.

4. After the HA container is empty, replace it with sodium chloride 0.9% and rinse the infusion system with up to 150 ml in order to recover the remaining albumin from the infusion system without changing the flow rate.

According to this use, up to 600 ml Albumin solution 20% can be filtered with up to 4 Hepalbin[®]-Adsorbents in a row, provided that the pre- and post-rinsing volumes are used according to Table 1.

4.3 Example C, for the use of 6 Hepalbin[®]-Adsorbents, e.g. for Single Pass Albumin Dialysis (SPAD)

This example describes the use of 6 Hepalbin[®]-Adsorbents for the purification of 600 to 800 ml Albumin solution 20% and is pictured in **Fig 1C**. In order to maintain the pressure for perfusion in an acceptable range, two columns of 3 Hepalbin[®]-Adsorbents in a row have to be connected by a “Y”-connector* for a parallel flow. The outlet of each “column” will be connected by one particle filter for each, which are converged using a 3-way stopcock as used for intravenous infusion systems.

Using an infusion pump without a drop counting control mechanism may be preferred, as albumin tends to cause foam in the dripping chamber which may cause false alarms due to the missing signal induced by the drops.

1. Connect part **1** (appr. 600 ml sodium chloride) with part **2** (infusion line with an infusion pump) and use an “Y”-connector (**6**) to enable the connection of two rows of 3 Hepalbin[®]-Adsorbents (**3**) each. Connect to each row a particle filter (**4**) and use a 3-way stopcock (**7**) to have the two particle filter lines converging as pictured in **Fig. 1C**.

2. During filling and pre-rinsing, the Hepalbin[®]-Adsorbent and the particle filter should be placed in a way that the fluid enters at the bottom and flows upwards. For achieving this, one could use a bottle holder (**Fig. 2**) and a clamp (**Fig. 3**) to fix the cluster temporarily on the infusion stand. This positioning of “upside down” ensures a complete and almost bubble free filling of the Hepalbin[®]-Adsorbents and the particle filters. The fluid leaving the particle filter should be discarded until 600 ml are used.

3. Replace the sodium chloride with the container of HA 20% solution. Connect the 3-way stopcock with the collection bag (**5**). Adjust the flow rate to approximately 900 ml/h. If the remaining sodium chloride in the system should not be used, from now, up to 350 ml can be discarded before the infusion line must be connected. This will reduce the dilution without risking significant loss of Albumin.

4. After proceeding the whole amount of HA solution, replace the container with sodium chloride 0.9% and rinse the infusion system with up to 350 ml in order to recover the remaining Albumin from the infusion system without changing the flow rate.

* Two “Y”-connectors are to be found in each packing unit at 20 Hepalbin[®]-Adsorbents.

4.4 Example D, for the use of 3 Hepalbin[®]-Adsorbents in a row to purify 80g of Albumin for Extracorporeal Albumin Dialysis (ECAD)

This example describes the use of 3 Hepalbin[®]-Adsorbents for the purification of 400 ml Albumin 20%. It is practically relevant in MARS[®] (Molecular Adsorbents Recirculating System). The use of 3 adsorbents is pictured in **Fig. 1D**.

Using an infusion pump without a drop counting control mechanism may be preferred, as albumin tends to cause foam in the dripping chamber which may cause false alarms due to the missing signal induced by the drops.

1. Connect part **1** (appr. 300 ml sodium chloride), part **2** (infusion line with an infusion pump) and 3 Hepalbin[®]-Adsorbents (**3**) in a row with the particle filter (**4**) as pictured in **Fig. 1D**.

2. During filling and pre-rinsing with the sodium chloride, the Hepalbin[®]-Adsorbent and the particle filter should be placed in a way that the fluid enters at the bottom and flows upwards. For achieving this, one could use a bottle holder (**Fig. 2**) and a clamp (**Fig. 3**) to fix the cluster temporarily on the infusion stand. This ensures a complete and almost bubble free filling of the Hepalbin[®]-Adsorbent and the particle filter with sodium chloride. The fluid leaving the particle filter should be discarded until about 300 ml are used, according to Table 1.

3. Replace the sodium chloride with 400 ml of HA 20% solution (If more than one container is needed, it can be replaced consecutively.). Adjust the flow rate to 450-600 ml/h. After about 100 ml HA are proceeded, connect the receptacle bag for recovering the albumin.

4. After the 400 ml HA container is empty, replace it with sodium chloride 0.9% and rinse the infusion system with 200 to 250 ml in order to recover the remaining albumin from the infusion system without changing the flow rate.

4.5 Technical Support Devices

The following technical Support Devices can be used for applying the Hepalbin®-Adsorbent:


Figure 2: infusion bottle holder

Figure 3: coated clamp with hook

- Infusion bottle holder [Fig. 2] to fix the cluster [Fig. 4] in particular (to be found in standard hospitals)
- clamp with hook [Fig. 3] to fix the Hepalbin®-Adsorbent to the infusion bottle holder [Fig. 4] or to hang a single adsorbent on the infusion stand (can be delivered by Albutec)
- „Y“-connector (delivered with a unit at 20) and 3-way stopcock (to be found in standard hospitals)

Figure 4: Hepalbin®-Adsorbent (cluster of six)

5 Limitations

In clinical applications, only approved pharmaceutical sodium chloride 0.9% and albumin solutions are allowed to be used in the application of the Hepalbin®-Adsorbent.

6 Warnings

The Hepalbin®-Adsorbent:

- is a single use device;
- should not be exposed to pressures higher than 1.2 bar;
- is pyrogen free and steril only if the primary package is intact and not damaged;
- can be used only if not expired;
- requires sterile handling.

7 Data

The Hepalbin®-Adsorbent is capable of reducing octanoate and acetyltryptophanate from commercial, pharmaceutical albumin preparations which contain those substances in inphysiological concentrations (5.3 mol/mol albumin!).

One Hepalbin®-Adsorbent is designed to reduce the octanoate/albumin (and correspondingly acetyltryptophanate/albumin) molecular ratio from 5.3 to less than 0.3. This is guaranteed by the use of 1 Hepalbin®-Adsorbent for 20 to 27 g albumin at flow rates of approximately 150 ml/h per Hepalbin®-Adsorbent if 20% solutions are used. The flow rate can be increased for saving time, but more than doubling the flow rate per adsorbent will reduce the reduction significantly. If using 3 Hepalbin®-Adsorbents and not more than 400 ml Albumin 20% solution for example, the optimal reduction will occur at flow rates of up to 3 times 150 ml/h equals 450 ml/h.


Figure 5: Octanoate concentration after proceeding 100 ml HA 20% per Hepalbin®-Adsorbent with different flow rates

8 Troubleshooting

One potential cause of poor performance is the insufficient filling of the device, which may cause air to remain in the Hepalbin®-Adsorbent and in the particle filter. This could cause

- a) ineffective use of the Hepalbin®-Adsorbent volume, which may cause less reduction of stabilizers
- b) ineffective use of the particle filter membrane area, which may cause higher pressure for the infusion pump or inability to infuse by gravity (if one adsorbent is used).

For prevention the adsorbent should be filled upside down that the fluid enters at the bottom. A small amount of bubbles in the Hepalbin®-Adsorbent is inevitable but not limiting. A slope of about 3 degree during filling is to the best advantage.

The air from the particle filter can be dispatched through the ventilation membrane by holding it up.

Under certain circumstances, air could enter the particle filter also via the "de-airing" access. This could be prevented by "clamping" the particle filter line when it is disconnected from the rinsing bags or the IV line.