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Please read this information carefully before use!

english

### General directions

Filter for 20% human albumin infusion solution (HA) and dialysate which removes undesired substances, such as stabilizers (octanoate/N-acetyl-tryptophanate), bile acids, bilirubin and creatinine.


**CAVE!** Infusion only with enclosed particle filter (RoweFil 48); suitable for use with infusion lines for infusion pumps.

**CAVE!** Use at least 100 ml of isotonic sodium chloride per Hepalbin®-Adsorbent for pre rinsing!

### Examples of application

For use of 1 filter in 100 ml 20% albumin no infusion pump is required. If more than 1 filter is used in a row, an infusion pump must be used.

As in example B you can put 2 up to 4 Hepalbin®-Adsorbents in a row. If using more Hepalbin®-Adsorbents you need to add an Y-connector (two enclosed in a packaging unit) in order to generate two columns (see example C).

Legend:

- 1 medical solution which shall be filtered
- 2 infusion lines 2a without, 2b with infusion pump
- 3 Hepalbin®-Adsorbent
- 4 particle filter RoweFil 48
- 5 A, B venous access; C receptacle bag for extracorporeal treatment
- 6 Y-connector
- 7 3-way stopcock

Table: Recommendations with respect to the quantity of HA

HA (20%) volume [ml]	100	200	300	400	500	600	700	800
Number of adsorbents	1	2	3	3	4	4	6	6
Flow rate [ml/h]	150	300	450	450	600	600	900	900
Minimum of sodium chloride for post-purge [ml]	100	150	200	200	250	250	350	350

### Preparation

- use a common metallic holder for infusion bottles for one or more Hepalbin®-Adsorbents
- assemble components as shown in figures (consider flow direction)
- avoiding bubbles, fill and rinse infusion system with sodium chloride (CAVE! Hepalbin®-Adsorbent and particle filter bottom up)
- during filling and rinsing, 100 ml sterile (infusion grade) 0.9% sodium chloride solution per adsorbent will enter the system out of which 30 ml to 50 ml will be discarded and 70 ml to 50 ml will remain in the system before connecting the albumin bottle
- discard excessive volume, attach HA, connect to patient

### Performance

- CAVE! Close infusion line below particle filter in case of interruption such as changing bottles
- 1 drop per second (to be monitored in drip chamber) is the optimal flow rate without using infusion pump
- the minimum of a rinsing volume to wash out albumin remains is given in the table
- during filling and rinsing, 100 ml sterile (infusion grade) 0.9% sodium chloride solution per adsorbent will enter the system out of which 30 ml will be discarded if only 1 Hepalbin®-Adsorbent is used, but if more adsorbents are used, an additional volume of 50 ml per adsorbent shall be discarded; the rest of the rinsing solution remains in the device
- for further directions see product information<sup>1)</sup>

### Disposal

Discard Hepalbin®-Adsorbent complying legal requirements for hazardous waste.

### Restrictions

Use only with Human Albumin and medical solutions which are approved as pharmaceuticals.

### Precautions

- for single use only (the ability of binding toxins is limited and not reversible, resterilisation results in material damage)!
- pressure-resistant up to 1,2 bar
- sterile and free of pyrogens only if package is intact
- do not use after expiration date
- handle with hygienic operating principles

### Performance characteristics

If applied as advised (table) there could be reached a ratio stabilizer/albumin well below 0.3 mol/mol.

### Sterilization

The Hepalbin®-Adsorbent is sterilized by gamma radiation.

### Storage

Storage of Hepalbin®-Adsorbent at room temperature, in dry conditions without direct sun exposure.