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
Please read information carefully before use !

english

### General notices

Filter cluster for the removal of undesired substances such as albumin stabilizers (octanoat, tryptophanat), bile acids, bilirubins and creatinine.


Application by trained staff only!

#### Application example A

Removal of undesired substances and toxins (i.e. octanoat, bile acids, bilirubins and creatinine) from toxin enriched dialysis and infusion solutions.

**CAVE!** Infusion solutions must be filtered after Hepalbin-Cluster<sup>12</sup> perfusion (pore size 0.22 µm).

#### Application example B

Removal of undesired substances and toxins (i.e. octanoat, bile acids, bilirubins and creatinine) from recirculating dialysis circuits containing hemofiltration solution or albumin dialysis solution.

*toxin permeable dialysator is required*

### Application notices

- Application example B requires the use of a toxin permeable dialysator (i.e. highflux dialysators are most likely suitable for the named substances).
- Application example B requires the use of a suitable dialysis platform.

### Warning notices

- The Hepalbin®-Cluster<sup>12</sup> is intended for single use only - the toxin adsorption capacity is limited and irreversible, resterilisation may cause material damage.
- The Hepalbin®-Cluster<sup>12</sup> is intended for pressures lower than 1 bar only.
- The Hepalbin®-Cluster<sup>12</sup> is sterile and free of pyrogens only if package is undamaged.
- The Hepalbin®-Cluster<sup>12</sup> is to be used before storage limit is reached (see identification label).
- The Hepalbin®-Cluster<sup>12</sup> requires hygienic and germ-free operation procedures.
- The Hepalbin®-Cluster<sup>12</sup> needs to be flushed by sterile saline solution of 1 liter minimum.

### Sterilisation

The Hepalbin®-Cluster<sup>12</sup> is sterilized via gamma radiation.

### Storage

The Hepalbin®-Cluster<sup>12</sup> needs to be stored by room temperature and in a dry environment out of direct sunlight.